

# آزمایشگاه کولایف تهران

## Control Chart

Date-Run	Value
1400/9/1-1	157
1400/9/3-1	154
1400/9/4-1	152
1400/9/6-1	151
1400/9/7-1	150.8
1400/9/8-1	147
1400/9/9-1	146
1400/9/10-1	147
1400/9/11-1	145
1400/9/13-1	144
1400/9/14-1	144
1400/9/15-1	153
1400/9/16-1	154
1400/9/20-1	147.4
1400/9/21-1	147.9
1400/9/22-1	155
1400/9/23-1	155
1400/9/24-1	151
1400/9/25-1	152

High					
mindray BS-800M	دستگاه	AST	تست	1400/9/1	تاریخ
1:2s	قوانین تفسیر	Control	Randox Hum Asy L3	Lot. No	1173UE

Mean	149	SD	8.5	Calculated :	Mean	150.16	SD	4.03	CV %	2.68 (6.15)	Sigma	5.94	N	19
------	-----	----	-----	--------------	------	--------	----	------	------	-------------	-------	------	---	----

