

آزمایشگاه پاتوبیولوژی و ژنتیک فرجاد

Control Chart

Date-Run	Value
1399/5/1-1	56
1399/5/2-1	55
1399/5/8-1	57
1399/5/8-2	53
1399/5/21-1	54
1399/5/21-2	56
1399/5/21-3	56
1399/5/21-4	60
1399/5/21-5	61
1399/5/21-6	60
1399/5/21-7	62
1399/5/27-1	61
1399/5/28-1	61
1399/6/1-2	59
1399/6/6-1	58
1399/6/11-1	60
1399/6/22-1	60
1399/6/25-1	55
1399/6/31-1	58
1399/7/8-3	59
1399/7/14-1	56
1399/7/15-1	57
1399/7/19-1	56
1399/8/1-3	59
1399/8/10-1	60
1399/8/11-2	57
1399/8/12-1	55
1399/8/15-2	56
1399/8/17-1	56
1399/8/18-1	58
1399/8/19-2	60
1399/8/20-1	57
1399/8/22-2	59
1399/8/24-1	61
1399/8/27-2	57

Normal					
mindray BS-800M		دستگاه	AST-ST	تست	1400/1/17 تاریخ
1:3s		قوانین تفسیر	Control	Delta D.P. Normal-SI	Lot. No 588

Mean	58.48	SD	1.70	Calculated :	Mean	58.36	SD	1.84	CV %	3.15	Sigma	6.19	N	116
------	-------	----	------	--------------	------	-------	----	------	------	------	-------	------	---	-----



Date-Run	Value
1399/8/28-2	56
1399/8/29-1	61
1399/9/1-4	60
1399/9/3-1	61
1399/9/4-1	58
1399/9/5-1	61
1399/9/6-2	58
1399/9/8-4	60
1399/9/9-1	61
1399/9/10-2	59
1399/9/11-3	58
1399/9/12-2	58
1399/9/15-3	61
1399/9/16-1	60
1399/9/17-2	59
1399/9/18-2	61
1399/9/19-3	58
1399/9/23-1	58
1399/9/24-2	58
1399/9/26-1	59
1399/9/27-3	59
1399/10/2-2	61
1399/10/8-5	60
1399/10/13-1	59
1399/10/21-8	61
1399/10/23-2	58
1399/10/24-3	58
1399/10/29-2	56
1399/10/30-1	58
1399/11/1-5	56
1399/11/4-2	56
1399/11/11-4	57
1399/11/12-3	56
1399/11/13-2	54
1399/12/4-3	58.4
1399/12/13-3	56
1399/12/16-2	58.7
1399/12/17-3	56
1399/12/18-3	60
1399/12/23-3	59
1399/12/25-2	58.8
1399/12/26-2	60
1399/12/27-4	60
1400/1/14-2	57.9
1400/1/15-4	57.4

Date-Run	Value
1400/1/16-3	56.8
1400/1/17-4	57.4
1400/1/18-4	56.1
1400/1/22-3	58.8
1400/1/28-3	57.2
1400/2/2-1	59.8
1400/2/6-1	60
1400/2/7-1	58.8
1400/2/7-1	58.8
1400/2/8-1	59.3
1400/2/9-1	61
1400/2/11-1	60
1400/2/13-1	58.9
1400/2/15-1	59
1400/2/18-1	59
1400/2/19-1	59.3
1400/2/20-1	57.3
1400/2/21-1	60
1400/2/22-1	58
1400/2/26-1	58.7
1400/2/27-2	57.8
1400/2/29-1	56.7
1400/3/3-1	59
1400/3/6-1	57.2
1400/3/8-1	58
1400/3/9-1	59.8
1400/3/10-1	57.1
1400/3/11-1	58.8
1400/3/12-1	58.8
1400/3/13-1	62
1400/3/19-1	59.2
1400/3/22-1	59.5
1400/3/24-1	58.9
1400/3/26-1	57.3
1400/3/29-1	57.4
1400/3/30-1	57