

آزمایشگاه کولایف تهران

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.69	
Mean	156.43	Mean	281.38
SD	2.36	SD	4.75
CV	1.51	CV	1.69
Day-Run	Value	Day-Run	Value
1 - 1	157	1 - 1	274
2 - 1	155	2 - 1	275
3 - 1	153	3 - 1	289
4 - 1	154	4 - 1	281
6 - 1	157	6 - 1	282
7 - 1	154	7 - 1	286
8 - 1	153	8 - 1	281
9 - 1	158	9 - 1	282
10 - 1	155	10 - 1	274
11 - 1	158	11 - 1	283
13 - 1	153	13 - 1	279
14 - 1	158	14 - 1	286
15 - 1	153	15 - 1	274
16 - 1	158	16 - 1	276
20 - 1	158	20 - 1	282
21 - 1	158	21 - 1	287
22 - 1	160	22 - 1	283
23 - 1	157	23 - 1	285
24 - 1	160	24 - 1	280
25 - 1	157	25 - 1	281
27 - 1	159	27 - 1	289

Westgard	تفسیر	CHOL	آزمایش	mindray BS-800M	دستگاه	9	ماه	1400	سال
Normal Control	Randox	Mean	157	SD	5.83	Exp. Date	1400/5/21	Lot. No	1489UN

