

برگشت (/IACLD/Panel/Labs/Answers/36743)

چاپ

دوره چهل و پنجم تابستان 1402

S2 Results

نمونه مجهول ToRCH

کد نمونه 1080	EQAP-19374	آزمایشگاه کولایف
---------------	------------	------------------

Anti-CMV (IgG)

Anti-HSV1 (Anti-Herpes Type1 IgG)

Anti-Rubella (IgG)

Anti-Toxoplasma (IgG)

Anti-CMV (IgG)

جدول شماره 1. توزیع فراوانی نتایج گزارش شده توسط شرکت کنندگان براساس روش

تذکر. نتیجه مورد قبول برای این آزمایش پاسخ + بوده است. نتیجه ارائه شده توسط آزمایشگاه شما با علامت ← مشخص شده است.

- ثبت نام کنندگان 222
- شرکت کنندگان 198
- کل درصد مشارکت 89.19 %
- پاسخهای دریافتی 177
- درصد پاسخدهی 89.39 %

Method	Count				Percent		
	Total	Positive	Indeterminate	Negative	Positive	Indeterminate	Negative
CL	57	55	0	2	%96.49	%0.00	%3.51
ECL	26	25	1	0	%96.15	%3.85	%0.00
ELFA	26	26	0	0	%100.00	%0.00	%0.00
ELISA	67	67	0	0	%100.00	%0.00	%0.00
Immunochromatography	0	0	0	0	%0.00	%0.00	%0.00
Other	1	1	0	0	%100.00	%0.00	%0.00

Method	Kit	Result	Count	Your Answer
CL	Abbott	Positive	8	
CL	Architect	Positive	6	
CL	Autobio Diagnostic	Positive	2	

Method	Kit	Result	Count	Your Answer
CL	Biomerieux	Positive	1	
CL	Chorus	Positive	1	
CL	Dia Sorin	Positive	24	
CL	Liaison	Not Reported	1	
CL	Liaison	Positive	4	
CL	Maglumi	Positive	3	
CL	Ortho	Positive	1	
CL	Roche	Positive	1	
CL	Siemens	Negative	2	
CL	Siemens	Positive	4	
ECL	Architect	Positive	2	
ECL	Biomerieux	Positive	1	
ECL	Cobas	Positive	2	
ECL	Dia Sorin	Positive	2	
ECL	ImmunoLab	Positive	1	
ECL	Liaison	Positive	1	
ECL	Maglumi	Positive	1	
ECL	Roche	Not Reported	2	
ECL	Roche	Indeterminate	1	
ECL	Roche	Positive	15	
ELFA	Biomerieux	Positive	21	
ELFA	Chorus	Positive	3	
ELFA	Diesse-Tec	Positive	1	
ELFA	Maglumi	Positive	1	
ELISA	Not Reported	Positive	3	
ELISA	Chorus	Positive	4	
ELISA	Diesse-Tec	Positive	1	
ELISA	Pishtaz Teb	Positive	52	←
ELISA	Siemens	Positive	1	
ELISA	Vircell	Positive	6	

Method	Kit	Result	Count	Your Answer
Other	Abbott	Positive	1	

Copyright © 2024 IACLD. All rights reserved.