

برگشت (/IACLD/Panel/Labs/Answers/36745)

چاپ

دوره چهل و هفتم زمستان 1402

S1.1 Results



نمونه مجهول ایمنوسرولوژی برای آزمایش HBsAg, HIVAb, S1.1 HCVAb

| | | |
|---------------|------------|------------------|
| کد نمونه 1826 | EQAP-19374 | آزمایشگاه کولایف |
|---------------|------------|------------------|

| | | |
|-------------|-------------|-------------|
| S1/1 HBs Ag | S1/1 HCV Ab | S1/1 HIV Ab |
|-------------|-------------|-------------|

S1/1 HBs Ag

جدول شماره 1. توزیع فراوانی نتایج گزارش شده توسط شرکت کنندگان براساس روش

تذکر. نتیجه مورد قبول برای این آزمایش پاسخ  بوده است. نتیجه ارائه شده توسط آزمایشگاه شما با علامت  مشخص شده است.

- ثبت نام کنندگان 1658
- شرکت کنندگان 1545
- کل درصد مشارکت 93.18 %
- پاسخهای دریافتی 1493
- درصد پاسخدهی 96.63 %

| Method | Count | | | | Percent | | |
|----------------------|-------|----------|---------------|----------|----------|---------------|----------|
| | Total | Positive | Indeterminate | Negative | Positive | Indeterminate | Negative |
| CL | 109 | 1 | 0 | 108 | %0.92 | %0.00 | %99.08 |
| ECL | 154 | 2 | 0 | 152 | %1.30 | %0.00 | %98.70 |
| ELFA | 58 | 0 | 0 | 58 | %0.00 | %0.00 | %100.00 |
| ELISA | 1143 | 97 | 11 | 1035 | %8.49 | %0.96 | %90.55 |
| Immunochromatography | 13 | 0 | 0 | 13 | %0.00 | %0.00 | %100.00 |
| Other | 13 | 1 | 0 | 12 | %7.69 | %0.00 | %92.31 |

| Method | Kit | Result | Count | Your Answer |
|--------|--------------------|----------|-------|-------------|
| CL | Not Reported | Negative | 1 | |
| CL | Abbott | Negative | 22 | |
| CL | Architect | Negative | 9 | |
| CL | Autobio Diagnostic | Negative | 1 | |
| CL | Dia Sorin | Negative | 21 | |

| Method | Kit | Result | Count | Your Answer |
|--------|--------------------------|--------------|-------|-------------|
| CL | Dia Sorin | Positive | 1 | |
| CL | ImmunoDiagnostic Systems | Negative | 1 | |
| CL | Liaison | Negative | 3 | |
| CL | Mindray | Negative | 19 | |
| CL | Ortho | Negative | 3 | |
| CL | Roche | Negative | 3 | |
| CL | Siemens | Negative | 25 | |
| ECL | Not Reported | Negative | 3 | |
| ECL | Abbott | Negative | 8 | |
| ECL | Architect | Negative | 3 | |
| ECL | Biomerieux | Negative | 2 | |
| ECL | Biomerieux | Positive | 1 | |
| ECL | Cobas | Negative | 12 | |
| ECL | Cobas | Positive | 1 | |
| ECL | Dia Sorin | Negative | 7 | |
| ECL | Electa Lab | Negative | 1 | |
| ECL | Liaison | Negative | 3 | |
| ECL | Mindray | Negative | 4 | |
| ECL | Roche | Negative | 95 | |
| ECL | Siemens | Negative | 14 | |
| ELFA | Not Reported | Negative | 3 | |
| ELFA | Biomerieux | Not Reported | 1 | |
| ELFA | Biomerieux | Negative | 49 | |
| ELFA | Cobas | Negative | 1 | |
| ELFA | Pishtaz Teb | Negative | 1 | |
| ELFA | Roche | Negative | 1 | |
| ELFA | Siemens | Negative | 1 | |
| ELFA | Tosoh | Negative | 2 | |
| ELISA | Not Reported | Not Reported | 1 | |
| ELISA | Not Reported | Negative | 36 | |
| ELISA | Not Reported | Positive | 8 | |
| ELISA | Abbott | Negative | 2 | |
| ELISA | Abbott | Positive | 1 | |

| Method | Kit | Result | Count | Your Answer |
|--------|-------------------------|---------------|-------|-------------|
| ELISA | Abon | Positive | 1 | |
| ELISA | Ai De | Positive | 1 | |
| ELISA | Apdia | Negative | 2 | |
| ELISA | Arya Mabna Tashkhis | Not Reported | 1 | |
| ELISA | Arya Mabna Tashkhis | Negative | 101 | |
| ELISA | Arya Mabna Tashkhis | Positive | 9 | |
| ELISA | Autobio Diagnostic | Positive | 1 | |
| ELISA | Binding Site | Negative | 1 | |
| ELISA | Biomerieux | Negative | 1 | |
| ELISA | Bio-Pharm | Negative | 1 | |
| ELISA | Dia Lab | Negative | 1 | |
| ELISA | Diagnostic | Negative | 2 | |
| ELISA | Diapro | Negative | 36 | |
| ELISA | Diapro | Positive | 1 | |
| ELISA | Gbc | Indeterminate | 3 | |
| ELISA | Gbc | Negative | 91 | |
| ELISA | Gbc | Positive | 4 | |
| ELISA | General Biologicals(GB) | Not Reported | 1 | |
| ELISA | General Biologicals(GB) | Indeterminate | 1 | |
| ELISA | General Biologicals(GB) | Negative | 32 | |
| ELISA | General Biologicals(GB) | Positive | 1 | |
| ELISA | Mindray | Negative | 3 | |
| ELISA | Monobind | Negative | 1 | |
| ELISA | Noyan Negin Parsian | Negative | 2 | |
| ELISA | Others | Negative | 2 | |
| ELISA | Padtan Teb | Negative | 2 | |
| ELISA | Pars Azmoon | Negative | 1 | |
| ELISA | Pishtaz Teb | Indeterminate | 7 | |
| ELISA | Pishtaz Teb | Negative | 689 | |
| ELISA | Pishtaz Teb | Positive | 68 | ← |
| ELISA | Wantai | Negative | 29 | |
| ELISA | Wantai | Positive | 2 | |

| Method | Kit | Result | Count | Your Answer |
|----------------------|----------------------------|--------------|-------|-------------|
| Immunochromatography | Not Reported | Negative | 1 | |
| Immunochromatography | Abon | Negative | 7 | |
| Immunochromatography | Addcare | Negative | 1 | |
| Immunochromatography | Arya Mabna Tashkhis | Negative | 1 | |
| Immunochromatography | Medical Biological Service | Negative | 1 | |
| Immunochromatography | Mindray | Negative | 1 | |
| Immunochromatography | Siemens | Negative | 1 | |
| Other | Not Reported | Negative | 3 | |
| Other | Abbott | Negative | 1 | |
| Other | Abon | Negative | 1 | |
| Other | Biomerieux | Not Reported | 1 | |
| Other | Biomerieux | Negative | 3 | |
| Other | Biomerieux | Positive | 1 | |
| Other | Diagnostic | Negative | 1 | |
| Other | Siemens | Negative | 2 | |
| Other | Tosoh | Negative | 1 | |
| RIA | Biomerieux | Negative | 1 | |
| RIA | Pishtaz Teb | Negative | 1 | |

Copyright © 2024 IACLD. All rights reserved.