

برگشت (/IACLD/Panel/Labs/Answers/36745)

چاپ

دوره چهل و هفتم زمستان 1402

B3 Results

نمونه مجهول پروتئین ادرار

کد نمونه 1795	EQAP-19374	آزمایشگاه کولایف
---------------	------------	------------------

عملکرد عالی (رنگ آبی)	DI	< 0.5
عملکرد مطلوب (رنگ سبز)	DI	0.5 - 1.0
نتیجه قابل قبول اما در حدود مرز (رنگ زرد)	DI	1.0 - 2.0
نیاز به بررسی روش و کنترل کالیبراسیون (رنگ نارنجی)	DI	2.0 - 3.0
وجود نقص و نیاز به پیگیری فوری (رنگ قرمز)	DI	> 3.0
نتیجه گزارش شده خارج از محدوده است (Out Of Range)	DI	> 10.0

برای مشاهده آنالیز آماری تمامی همگروه ها، در جدول زیر بر روی نام آنالیت مورد نظر کلیک کنید.

Analyte	Unit	Lab Result	PeerGroup Count	PeerGroup Mean	PeerGroup %CV	Lab DI
Urine Protein (/IACLD/Panel/Analysis/Report/36745/595/183771)	mg/L	350.00	237	368.73	9.67	0.51
Protein/Cr Ratio (/IACLD/Panel/Analysis/Report/36745/595/183802)	mg/g Cr	581.40				
Protein (24h) (/IACLD/Panel/Analysis/Report/36745/595/183824)	mg/24h	700.00				
Creatinine B3 (/IACLD/Panel/Analysis/Report/36745/595/183833)	mg/dL	60.20	358	63.97	12.22	0.59
Creatinine 24 (/IACLD/Panel/Analysis/Report/36745/595/183868)	mg/24h	1204.00				
Calcium (/IACLD/Panel/Analysis/Report/36745/595/183876)	mg/dL	4.30	406	4.52	10.19	0.48
Calcium 24 (/IACLD/Panel/Analysis/Report/36745/595/183908)	mg/24h	86.00				
Calcium/Creatinine Ratio (/IACLD/Panel/Analysis/Report/36745/595/183927)	mg/mg	0.07				
Uric Acid B3 (/IACLD/Panel/Analysis/Report/36745/595/183938)	mg/dL	16.40	58	15.74	20.53	0.42
Uric Acid 24 (/IACLD/Panel/Analysis/Report/36745/595/183962)	mg/24h	328.00				

UricAcid/Creatinine Ratio (/IACLD/Panel/Analysis/Report/36745/595/183978)	mg/mg	0.27				
---	-------	------	--	--	--	--

Copyright © 2024 IACLD. All rights reserved.