

برگشت (/IACLD/Panel/Labs/Answers/36744)

چاپ

دوره چهل و ششم پائیز 1402

## S1.1 Results

## نمونه مجهول ایمنوسرولوژی برای آزمایش HBsAg, HIVAb, S1.1 HCVAb

کد نمونه 1805	EQAP-19374	آزمایشگاه کولایف
---------------	------------	------------------

S1/1 HBs Ag

S1/1 HCV Ab

S1/1 HIV Ab

## S1/1 HBs Ag

جدول شماره 1. توزیع فراوانی نتایج گزارش شده توسط شرکت کنندگان براساس روش

تذکر. نتیجه مورد قبول برای این آزمایش پاسخ + بوده است. نتیجه ارائه شده توسط آزمایشگاه شما با علامت ← مشخص شده است.

- ثبت نام کنندگان 1638
- شرکت کنندگان 1558
- کل درصد مشارکت 95.12 %
- پاسخهای دریافتی 1515
- درصد پاسخدهی 97.24 %

Method	Count				Percent		
	Total	Positive	Indeterminate	Negative	Positive	Indeterminate	Negative
CL	124	123	0	1	%99.19	%0.00	%0.81
ECL	142	139	0	3	%97.89	%0.00	%2.11
ELFA	76	74	0	2	%97.37	%0.00	%2.63
ELISA	1145	1081	15	49	%94.41	%1.31	%4.28
Immunochromatography	11	5	1	5	%45.45	%9.09	%45.45
Other	12	10	0	2	%83.33	%0.00	%16.67

Method	Kit	Result	Count	Your Answer
CL	Not Reported	Positive	3	
CL	Abbott	Positive	25	
CL	Architect	Positive	10	

Method	Kit	Result	Count	Your Answer
CL	Biomerieux	Positive	2	
CL	Cobas	Positive	1	
CL	Dia Sorin	Positive	22	
CL	ImmunoLab	Positive	1	
CL	Liaison	Positive	7	
CL	Mindray	Positive	17	
CL	Ortho	Positive	3	
CL	Siemens	Negative	1	
CL	Siemens	Positive	31	
CL	Vircell	Positive	1	
ECL	Not Reported	Positive	4	
ECL	Abbott	Positive	5	
ECL	Architect	Positive	1	
ECL	Biomerieux	Positive	4	
ECL	Cobas	Negative	1	
ECL	Cobas	Positive	9	
ECL	Dia Sorin	Positive	4	
ECL	Liaison	Positive	5	
ECL	Mindray	Positive	3	
ECL	Roche	Not Reported	1	
ECL	Roche	Negative	1	
ECL	Roche	Positive	95	
ECL	Siemens	Negative	1	
ECL	Siemens	Positive	8	
ECL	Tosoh	Positive	1	
ELFA	Anti Human	Positive	1	
ELFA	Biomerieux	Negative	2	
ELFA	Biomerieux	Positive	66	
ELFA	Others	Positive	1	
ELFA	Sam Fan Avaran	Positive	1	

Method	Kit	Result	Count	Your Answer
ELFA	Siemens	Positive	2	
ELFA	Tosoh	Positive	3	
ELISA	Not Reported	Not Reported	1	
ELISA	Not Reported	Indeterminate	3	
ELISA	Not Reported	Negative	2	
ELISA	Not Reported	Positive	36	
ELISA	Ai De	Positive	2	
ELISA	Anti Human	Positive	1	
ELISA	Arya Mabna Tashkhis	Not Reported	1	
ELISA	Arya Mabna Tashkhis	Indeterminate	4	
ELISA	Arya Mabna Tashkhis	Negative	3	
ELISA	Arya Mabna Tashkhis	Positive	78	
ELISA	Autobio Diagnostic	Positive	1	
ELISA	Bio Kit	Positive	1	
ELISA	Biomerieux	Positive	2	
ELISA	Bio-Pharm	Positive	2	
ELISA	Diapro	Negative	2	
ELISA	Diapro	Positive	31	
ELISA	Gbc	Negative	6	
ELISA	Gbc	Positive	67	
ELISA	General Biologicals(GB)	Negative	3	
ELISA	General Biologicals(GB)	Positive	31	
ELISA	Monobind	Positive	2	
ELISA	Padtan Teb	Negative	1	
ELISA	Padtan Teb	Positive	1	
ELISA	Pishtaz Teb	Not Reported	3	
ELISA	Pishtaz Teb	Indeterminate	7	
ELISA	Pishtaz Teb	Negative	28	
ELISA	Pishtaz Teb	Positive	735	
ELISA	Surase B-96(Tmb)	Positive	2	

Method	Kit	Result	Count	Your Answer
ELISA	Vironostika	Positive	1	
ELISA	Wantai	Indeterminate	1	
ELISA	Wantai	Negative	4	
ELISA	Wantai	Positive	88	←
Immunochromatography	Not Reported	Negative	1	
Immunochromatography	Not Reported	Positive	1	
Immunochromatography	Abon	Negative	4	
Immunochromatography	Abon	Positive	1	
Immunochromatography	Arya Mabna Tashkhis	Positive	1	
Immunochromatography	Biomerieux	Positive	1	
Immunochromatography	Medical Biological Service	Indeterminate	1	
Immunochromatography	Roche	Positive	1	
Other	Not Reported	Negative	1	
Other	Not Reported	Positive	1	
Other	Abbott	Positive	1	
Other	Acon	Negative	1	
Other	Biomerieux	Positive	4	
Other	Ortho	Positive	1	
Other	Others	Positive	1	
Other	Siemens	Positive	2	
RIA	Biomerieux	Positive	2	
RIA	Wantai	Positive	1	

Copyright © 2024 IACLD. All rights reserved.