

آزمایشگاه کولایف تهران

Control Chart

Date-Run	Value
1400/8/1-1	40.44
1400/8/3-1	39.89
1400/8/4-1	40.8
1400/8/5-1	40.36
1400/8/6-1	40.61
1400/8/8-1	38.71
1400/8/9-1	39.03
1400/8/10-1	39.46
1400/8/11-1	36.65
1400/8/13-1	41.81
1400/8/16-1	42.17
1400/8/17-1	40.46
1400/8/18-1	37.07
1400/8/19-1	42.79
1400/8/20-1	41.42
1400/8/22-1	39.61
1400/8/23-1	40.55
1400/8/25-1	42.16
1400/8/26-1	40.32
1400/8/27-1	40.74

High					
Architect	دستگاه	FSH	تست	1400/8/1	تاریخ
1:3s	قوانین تفسیر	Control	multichem	Lot. No	33402200

Mean	39.66	SD	1.55	Calculated :	Mean	40.25	SD	1.57	CV %	3.89 (5.5)	Sigma	6.40	N	20
------	-------	----	------	--------------	------	-------	----	------	------	------------	-------	------	---	----

