

آزمایشگاه کولایف تهران

Control Chart

Date-Run	Value
1400/8/1-1	14.57
1400/8/3-1	15.17
1400/8/4-1	15.24
1400/8/5-1	15.1
1400/8/6-1	15.21
1400/8/8-1	14.8
1400/8/9-1	14.93
1400/8/10-1	14.72
1400/8/11-1	14.68
1400/8/13-1	14.71
1400/8/15-1	15
1400/8/16-1	15.37
1400/8/17-1	15.16
1400/8/18-1	14.72
1400/8/19-1	15.08
1400/8/20-1	14.98
1400/8/22-1	14.56
1400/8/23-1	15.63
1400/8/24-1	15.47
1400/8/25-1	14.88
1400/8/26-1	15.32
1400/8/27-1	15.61

Normal					
Architect	دستگاه	FSH	تست	1400/8/1	تاریخ
1:3s - R:4s - 2of3:2s	قوانین تفسیر	Control	multichem	Lot. No	33402200

Mean	14.88	SD	0.64	Calculated :	Mean	15.04	SD	0.32	CV %	2.13 (5.5)	Sigma	11.67	N	22
------	-------	----	------	--------------	------	-------	----	------	------	------------	-------	-------	---	----

