

# آزمایشگاه کولایف تهران

## Control Serum - Levey Jennings Chart

Normal		High - CV : 4.58	
Mean	18.09	Mean	25.41
SD	0.68	SD	1.16
CV	3.77	CV	4.58
Day-Run	Value	Day-Run	Value
27 - 1	18.29	27 - 1	26.85
29 - 1	18.22	29 - 1	25.97
30 - 1	18.78	30 - 1	27.15
31 - 1	19.25	31 - 1	26.19
1 - 1	18.98	1 - 1	27.32
3 - 1	18.94	3 - 1	25.36
5 - 1	17.87	5 - 1	25.89
7 - 1	17.11	7 - 1	24.8
5 - 9	18.43	10 - 1	26
5 - 12	17.6	5 - 12	24.51
5 - 13	17.64	5 - 13	25.43
14 - 1	17.53	5 - 14	24.57
16 - 1	18.15	16 - 1	23.18
17 - 1	16.81	17 - 1	23.91
18 - 1	17.89	18 - 1	24.71
19 - 1	17.94	19 - 1	24.71

WHO	تفسیر	Homocysteine	آزمایش	Architect	دستگاه	6	ماه	1400	سال
Normal Control	Siemens	Mean	17.3	SD	1.73	Exp. Date	1400/3/17	Lot. No	33402200

