

آزمایشگاه کولایف تهران

Control Serum - Levey Jennings Chart

Normal	
Mean	5.86
SD	0.07
CV	1.2
Day-Run	Value
4 - 1	5.9
6 - 1	5.7
7 - 1	5.9
8 - 1	5.84
9 - 1	5.9
10 - 1	5.9
11 - 1	5.89
13 - 1	5.82
14 - 1	5.9
15 - 1	5.79
16 - 1	5.8
20 - 1	5.88
21 - 1	6
22 - 1	5.8
23 - 1	5.9
24 - 1	5.9
25 - 1	5.83
27 - 1	5.75

High - CV : 1.16	
Mean	4.7
SD	0.05
CV	1.16
Day-Run	Value
4 - 1	4.68
6 - 1	4.6
7 - 1	4.8
8 - 1	4.7
9 - 1	4.7
10 - 1	4.69
11 - 1	4.63
13 - 1	4.69
14 - 1	4.73
15 - 1	4.72
16 - 1	4.7
20 - 1	4.73
21 - 1	4.8
22 - 1	4.7
23 - 1	4.7
25 - 1	4.69
27 - 1	4.6

Westgard	تفسیر	t-pro	آزمایش	mindray BS-800M	دستگاه	9	ماه	1400	سال
Normal Control	Randox	Mean	5.83	SD	0.15	Exp. Date	1401/1/16	Lot. No	1489UN

