

آزمایشگاه کولایف تهران

Control Chart

Date-Run	Value
1400/9/1-1	277
1400/9/3-1	283
1400/9/4-1	287
1400/9/6-1	279
1400/9/8-1	280
1400/9/9-1	282
1400/9/10-1	280
1400/9/11-1	284
1400/9/13-1	278
1400/9/14-1	279
1400/9/15-1	276
1400/9/16-1	276
1400/9/20-1	283
1400/9/21-1	287
1400/9/22-1	282
1400/9/23-1	281
1400/9/24-1	280
1400/9/25-1	277
1400/9/27-1	278

High					
mindray BS-800M	دستگاه	TG	تست	1400/9/1	تاریخ
1:3s	قوانین تفسیر	Control	Randox Hum Asy L3	Lot. No	1173UE

Mean	280	SD	9.8	Calculated :	Mean	280.47	SD	3.31	CV %	1.18 (9.95)	Sigma	11.45	N	19
------	-----	----	-----	--------------	------	--------	----	------	------	-------------	-------	-------	---	----

