

# آزمایشگاه کولایف تهران

## Control Serum - Levey Jennings Chart

Normal		High		Low	
Mean	8.77	Mean	17.08	Mean	4.97
SD	0.73	SD	1.4	SD	0.52
CV	8.33	CV	8.18	CV	10.37
Day-Run	Value	Day-Run	Value	Day-Run	Value
3 - 1	9.4	3 - 1	16.3	3 - 1	5
4 - 1	8.7	4 - 1	16.6	4 - 1	5
6 - 1	10.7	6 - 1	15.3	6 - 1	6.1
7 - 1	8	7 - 1	15.9	7 - 1	4.6
8 - 1	8.6	8 - 1	18.6	8 - 1	5.1
9 - 1	9.4	9 - 1	19.6	9 - 1	5.9
10 - 1	8.3	10 - 1	15.7	10 - 1	4.9
11 - 1	8.6	11 - 1	19.6	11 - 1	4.9
13 - 1	7.9	13 - 1	17	13 - 1	4.8
14 - 1	8.3	14 - 1	17.1	14 - 1	4.9
15 - 1	8.6	15 - 1	16.7	15 - 1	4.4
16 - 1	8.7	16 - 1	16.2	16 - 1	4.7
17 - 1	8.8	17 - 1	17.5	17 - 1	4.3

Westgard	تفسیر	Tacrolimus	آزمایش	Cobas Immunoanalyzer	دستگاه	6	ماه	1400	سال
Normal Control	WBT	Mean	8.52	SD	0.85	Exp. Date	1400/5/25	Lot. No	
High Control	WBT	Mean	18	SD	1.8	Exp. Date	1400/5/25	Lot. No	
Low Control	WBT	Mean	5.01	SD	0.5	Exp. Date	1400/5/25	Lot. No	

