

آزمایشگاه کولایف تهران

Control Serum - Levey Jennings Chart

Normal		High		Low	
Mean	15.06	Mean	28.89	Mean	0.76
SD	0.37	SD	0.83	SD	0.04
CV	2.44	CV	2.87	CV	5.24
Day-Run	Value	Day-Run	Value	Day-Run	Value
20 - 1	14.47	20 - 1	27.65	20 - 1	0.71
21 - 1	15.07	21 - 1	29.83	21 - 1	0.73
23 - 1	14.94	23 - 1	28.85	23 - 1	0.8
24 - 1	15.61	24 - 1	27.88	24 - 1	0.77
25 - 1	15.3	25 - 1	28.85	25 - 1	0.74
26 - 1	15.16	26 - 1	27.98	26 - 1	0.76
30 - 1	15.48	30 - 1	29.69	30 - 1	0.75
31 - 1	14.91	31 - 1	28.98	31 - 1	0.8
1 - 1	14.53	1 - 1	29.24	1 - 1	0.84
2 - 1	15.14	2 - 1	29.97	2 - 1	0.73

Westgard	تفسیر	Testosterone	آزمایش	Architect	دستگاه	6	ماه	1400	سال
Normal Control	multichem	Mean	15.17	SD	0.63	Exp. Date	1400/1/25	Lot. No	33402200
High Control	multichem	Mean	28.5	SD	0.98	Exp. Date	1400/1/25	Lot. No	33402200
Low Control	multichem	Mean	0.74	SD	0.06	Exp. Date	1400/1/25	Lot. No	33402200

