

# آزمایشگاه پاتوبیولوژی و ژنتیک فرجاد

## Control Serum - Levey Jennings Chart

Normal		High		Low	
Mean	9.33	Mean	15.41	Mean	1.92
SD	0.39	SD	0.97	SD	0.22
CV	4.21	CV	6.29	CV	11.61
Day-Run	Value	Day-Run	Value	Day-Run	Value
1 - 1	10	1 - 1	13.6	3 - 1	2
2 - 1	10.1	2 - 1	13.9	8 - 1	1.7
12 - 1	8.9	13 - 1	15.6	12 - 1	2
13 - 1	9	18 - 1	15.3	22 - 1	1.7
18 - 1	9.3	19 - 1	14.7	27 - 1	2.22
19 - 1	9.3	20 - 1	14.9		
20 - 1	9.1	23 - 1	15.4		
22 - 1	9.1	24 - 1	15.6		
23 - 1	9.2	26 - 1	16.4		
24 - 1	9.5	27 - 1	16.7		
26 - 1	9.4	29 - 1	16.3		
30 - 1	8.8	30 - 1	15.3		
31 - 1	9.6	31 - 1	16.6		

Westgard	تفسیر	T4		آزمایش	ADVIA Immunoassay		دستگاه	3	ماه	1400	سال
Normal Control	Indox Premium Plus	Mean	9.16	SD	0.63	Exp. Date	1400/2/30	Lot. No	1765EC		
High Control	Indox Premium Plus	Mean	14.9	SD	1.08	Exp. Date	1400/2/30	Lot. No	1130EC		
Low Control	Indox Premium Plus	Mean	1.76	SD	0.21	Exp. Date	1400/2/30	Lot. No	1752EC		

