

آزمایشگاه کولایف تهران

Control Chart

Date-Run	Value
1400/9/1-1	126.5
1400/9/3-1	134
1400/9/6-1	128.4
1400/9/7-1	130.4
1400/9/8-1	126
1400/9/9-1	131
1400/9/10-1	127
1400/9/11-1	131
1400/9/13-1	126
1400/9/14-1	128.8
1400/9/15-1	126
1400/9/16-1	131
1400/9/20-1	126
1400/9/21-1	127.4
1400/9/22-1	130
1400/9/23-1	126.5
1400/9/24-1	127
1400/9/25-1	126
1400/9/27-1	127

High					
mindray BS-800M	دستگاه	Urea	تست	1400/9/1	تاریخ
1:3s	قوانین تفسیر	Control	Randox Hum Asy L3	Lot. No	1173UE

Mean	128	SD	5.15	Calculated :	Mean	128.21	SD	2.36	CV %	1.84 (6.05)	Sigma	8.11	N	19
------	-----	----	------	--------------	------	--------	----	------	------	-------------	-------	------	---	----

