

آزمایشگاه کولایف تهران

Control Chart

Date-Run	Value
1400/9/1-1	45
1400/9/3-1	44.2
1400/9/6-1	45.3
1400/9/7-1	45
1400/9/8-1	43.23
1400/9/9-1	46
1400/9/10-1	43.6
1400/9/11-1	44.47
1400/9/13-1	44.8
1400/9/14-1	43
1400/9/15-1	43
1400/9/16-1	45
1400/9/20-1	45
1400/9/21-1	44.6
1400/9/22-1	43.9
1400/9/23-1	44.9
1400/9/24-1	45
1400/9/25-1	43
1400/9/27-1	44

Normal					
mindray BS-800M	دستگاه	Urea	تست	1400/9/1	تاریخ
1:3s - R:4s - 2of3:2s	قوانین تفسیر	Control	Randox Hum Asy L2	Lot. No	1489UN

Mean	44.23	SD	1.9	Calculated :	Mean	44.37	SD	0.88	CV %	1.98 (6.05)	Sigma	7.58	N	19
------	-------	----	-----	--------------	------	-------	----	------	------	-------------	-------	------	---	----

