

آزمایشگاه کولایف تهران

Control Serum - Levey Jennings Chart

Normal	
Mean	81.83
SD	3.64
CV	4.45
Day-Run	Value
2 - 1	82.7
9 - 2	77
9 - 4	75.8
9 - 6	83.3
9 - 7	84.4
9 - 8	76.3
9 - 9	82.8
9 - 10	80.4
9 - 11	85.7
9 - 13	82.3
9 - 15	81.4
9 - 16	83.8
9 - 17	87.9

High - CV : 3.86	
Mean	154.61
SD	5.96
CV	3.86
Day-Run	Value
2 - 1	150.8
9 - 2	146.3
9 - 4	149.1
9 - 6	156.8
9 - 7	163.1
9 - 8	142
9 - 9	154.7
9 - 10	148.1
9 - 11	159
9 - 13	159.8
9 - 15	147.7
9 - 16	160.6
9 - 18	156.1
9 - 20	160.3
9 - 21	156.1
9 - 23	150.7
9 - 24	158.2
9 - 25	156.6
9 - 27	161.6

Westgard	تفسیر	Vitamin A	آزمایش	Agilent UV	دستگاه	9	ماه	1400	سال
Normal Control	RECIPE	Mean	83.73	SD	5.1	Exp. Date	1401/11/12	Lot. No	8870-8873

