

# آزمایشگاه کولایف تهران

## Control Serum - Levey Jennings Chart

High	
Mean	77.26
SD	1.88
CV	2.44
Day-Run	Value
9 - 2	78.5
9 - 4	77.5
9 - 7	80.8
9 - 13	77.8
9 - 15	77.2
9 - 17	76.8
9 - 20	74
9 - 22	77.2
9 - 24	75.5

Low	
Mean	20.89
SD	0.7
CV	3.36
Day-Run	Value
9 - 1	21.6
9 - 3	21.7
9 - 6	20.1
9 - 8	21.5
9 - 9	21.6
9 - 10	20.5
9 - 11	20.5
9 - 14	21.5
9 - 16	21.1
9 - 18	19.7
9 - 21	20.5
9 - 23	20
9 - 25	21.3

WHO	تفسیر	Vit D3	آزمایش	HPLC- HP	دستگاه	9	ماه	1400	سال
High Control	hChec-Recip AminA	Mean	77.76	SD	2.66	Exp. Date	1400/1/18	Lot. No	1289
Low Control	hChec-Recip AminA	Mean	20.76	SD	0.84	Exp. Date	1400/1/18	Lot. No	1289

